


## Synthesis of *tert*-butyl chloride – Protocol S




**Reaction.** In an Erlenmeyer flask, cool 17 mL (204 mmol) of concentrated hydrochloric acid (about 278% excess) with an ice-water bath to close to 0 °C. Transfer the cold acid to a separatory funnel and carefully add 4 g (54 mmol) of *tert*-butyl alcohol. Shake the mixture well, carefully venting. Swirl it occasionally for 20 minutes.

**Isolation.** Allow the layers to cleanly separate and remove the aqueous acid layer. Wash the product layer first with 10 mL of water, and then with 10 mL of 5% sodium hydrogen carbonate solution. Finally, wash the product layer with 10 mL of saturated sodium chloride solution, always removing the aqueous layer. Dry the product over anhydrous sodium sulphate for at least 10 minutes. Decant the product.

**Purification.** Purify the product using distillation.

**Safety.** See hazards associated with the reagents in Table 1.

**Greenness Assessment.** The evaluation was performed using the Green Star (GS) and the results are shown in Figure 1.



**Figure 1.** Greenness assessment (GS) for the synthesis of *tert*-butyl chloride

## Construction of the GS




Table 1 presents the hazards and scores associated with the substances involved and Table 2 presents the scores used to construct the green stars.

**Table 1.** Hazards for the synthesis of *tert*-butyl chloride, protocol S<sup>α</sup>

Substances involved	Step			Hazard code	Score: hazards to...		
	R	I	Pu		HH	E	P
<b>Stoichiometric reagents</b>							
Hydrochloric acid (CAS 7647-01-0)	✓			H314, H335	3	1	1
<i>tert</i> -Butyl alcohol (CAS 75-65-0)	✓			H225, H319, H332, H335	2	1	3
<b>Auxiliary substances</b>							
<b>Solvents</b>							
Sodium chloride (saturated solution)		✓		-	1	1	1
Sodium hydrogen carbonate (dilute solution)		✓		-	1	1	1
Water <sup>a,b</sup>		✓		-	1	1	1
<b>Other auxiliary substances</b>							
Anhydrous sodium sulphate (CAS 7757-82-6)		✓		-	1	1	1
<b>Product</b>							
<i>tert</i> -Butyl chloride (507-20-0)	✓	✓	✓	H225	1	1	3
<b>Waste</b>							
Carbon dioxide		✓		H280	1	1	2
Hydrochloric acid (dilute solution)		✓		-	1	1	1
Isobutylene	✓			H220, H280	1	1	3
Sodium chloride (aqueous solution)		✓		-	1	1	1
Sodium sulphate		✓		-	1	1	1
Water <sup>a,b</sup>		✓	✓	-	1	1	1

<sup>α</sup> R – Reaction; I – Isolation; Pu – Purification; HH – Human Health; E – Environment; P – Physical

<sup>a</sup> Renewable; <sup>b</sup> Degradable to innocuous products

**Table 2.** Scores used to construct the green star for the synthesis of *tert*-butyl chloride, protocol S<sup>a</sup>

Green Chemistry Principle	Reaction		Isolation		Purification		Global	
	s	Explanation	s	Explanation	s	Explanation	s	Explanation
<b>P1</b> Prevention	3	Isobutylene	3	Waste is innocuous	3	Waste is innocuous	3	Waste is innocuous
<b>P2</b> Atom Economy	1	Excess of hydrochloric acid > 10%, formation of by-products		NA		NA	1	Excess of hydrochloric acid > 10%, formation of by-products
<b>P3</b> Less hazardous chemical synthesis	1	Hydrochloric acid, H314		NA		NA	1	Hydrochloric acid, H314
<b>P5</b> Safer solvents and auxiliary substances	3	Solvents and auxiliary substances are not used	3	Solvents and auxiliary substances are innocuous	3	Solvents and auxiliary substances are not used	3	Solvents and auxiliary substances are innocuous
<b>P6</b> Increase energy efficiency	2	0 °C ≤ T ≤ 100 °C	3	Room temperature	2	0 °C ≤ T ≤ 100 °C	2	0 °C ≤ T ≤ 100 °C
<b>P7</b> Use renewable feedstocks	1	Substances not renewable	1	Substances not renewable	1	Substances not renewable	1	Substances not renewable
<b>P8</b> Reduce derivatives	3	One stage		NA		NA	2	Two stages
<b>P9</b> Catalysts	3	Without catalysts		NA		NA	3	Without catalysts
<b>P10</b> Design for degradation	1	Substances not degradable	1	Substances not degradable	1	Substances not degradable	1	Substances not degradable
<b>P12</b> Safer chemistry for accident prevention	1	Hydrochloric acid, H314, <i>tert</i> -butyl alcohol and <i>tert</i> -butyl chloride, H225, isobutylene, H220	1	<i>tert</i> -Butyl chloride, H225	1	<i>tert</i> -Butyl chloride, H225	1	Hydrochloric acid, H314, <i>tert</i> -butyl alcohol and <i>tert</i> -butyl chloride, H225, isobutylene, H220

<sup>a</sup>s – Score; NA – Not applicable

## References

Georgia Regents University,  
[http://www.aug.edu/~tcrute/lab\\_experiments/3411%20lab%20SN1%20prep%20of%20tertbutyl%20chloride.doc](http://www.aug.edu/~tcrute/lab_experiments/3411%20lab%20SN1%20prep%20of%20tertbutyl%20chloride.doc) (accessed April 2013).